


DATA
OUTPUT

ELECTRONIC BALANCES



8301-3000

- Automatic zero when power on
- Unit: g (gram), ct (carat), ozt (troy ounce)
- Weighing, counting, taring
- Overload alarm
- Data output



Code	8301-300	8301-600	8301-1200	8301-3000	8301-6000		
Maximum weighing	300g	600g	1200g	3000g	6000g		
Minimum weighing	0.2g	0.2g	0.2g	2g	2g		
Resolution	0.01g	0.01g	0.01g	0.1g	0.1g		
Weighing pan size	Ø135mm			190×160mm			
Verification scale interval	10d						
Working temperature	0~40℃						
Working humidity	≤80% RH						
Power supply	power adapter						
Dimension (L×W×H)	215×196×62mm						



calibration weight
(included in 8301-300,
8301-600 and 8301-1200)



windshield
(included in 8301-300,
8301-600 and 8301-1200)

OPTIONAL ACCESSORY

RS232 CABLE	8301-CABLE
-------------	------------




WEIGHING SCALES (OIML CERTIFICATE)



8002-30M

- LCD display, with backlight
- Zero-setting, subtracting tare value
- Overload alarm
- Automatic power off
- OIML certificate (No. 0200-NAWI-04601)



SPECIFICATION


Code	8002-6M	8002-15M	8002-30M	8002-6MD	8002-15MD	8002-30MD			
Maximum capacity	3/6kg	6/15kg	15/30kg	3/6kg	6/15kg	15/30kg			
Minimum capacity	20g	40g	100g	20g	40g	100g			
Division	1/2g	2/5g	5/10g	1/2g	2/5g	5/10g			
Verification interval	2g	5g	10g	2g	5g	10g			
Output	no			RS232					
Operating temperature	-10~40℃								
Power	rechargeable battery (6V/3.2Ah), power adapter (12V/500mA)								
Dimensions of pan	263×204mm								

WEIGHING SCALES (HIGH PRECISION)

- LCD display, with backlight
- Zero-setting, subtracting tare value
- Overload alarm
- Automatic power off



8001-30



SPECIFICATION

Code	8001-6	8001-15	8001-30	8001-6D	8001-15D	8001-30D
Maximum capacity	6kg	15kg	30kg	6kg	15kg	30kg
Minimum capacity	4g	10g	20g	4g	10g	20g
Division	0.2g	0.5g	1g	0.2g	0.5g	1g
Verification interval	2g	5g	10g	2g	5g	10g
Units	kg, lb			kg		
Output	no			RS232		
Operating temperature	-10~40℃					
Power	rechargeable battery (6V/3.2Ah), power adapter (12V/500mA)					
Dimensions of pan	263×204mm					



printer (optional)

OPTIONAL ACCESSORY (for scales with RS232)


Printer	8001-PRINTER
---------	--------------

COUNTING SCALES

- LCD display, with backlight
- Zero-setting, subtracting tare value, counting, accumulation, checking alarm
- Overload alarm
- Units: kg, g, lb, tj, hj, lb-oz



8101-30



SPECIFICATION

Code	8101-6	8101-15	8101-30	8101-6D	8101-15D	8101-30D			
Maximum capacity	6kg	15kg	30kg	6kg	15kg	30kg			
Minimum capacity	4g	10g	20g	4g	10g	20g			
Division	0.2g	0.5g	1g	0.2g	0.5g	1g			
Verification interval	2g	5g	10g	2g	5g	10g			
Output	no			RS232					
Operating temperature	-10~40℃								
Power	rechargeable battery (6V/3.2Ah), power adapter (12V/500mA)								
Dimensions of pan	263×204mm								



printer (optional)

OPTIONAL ACCESSORY (for scales with RS232)

Printer	8001-PRINTER
---------	--------------